


# BRIO DYNAMIC ULTRASOUND

## With an innovative Dynamic Mode in Ultrasound system, Brio touches new horizon of technology.

Dynamic mode uses penetration depth advantage of 1 & 3 Mhz along with power modulation to give layer by layer treatment.



## BRIO INTRODUCES HANDS FREE ULTRASOUND APPLICATOR SYSTEM :

- It will provide Precise and Targeted treatment
- Has special modulation program to avoid Hot Spots
- No need to move applicator, which saves operator's time & energy
- Zero operator errors for more efficient Treatment



Manufactured using world's best German piezo crystal

Ultrasound watch Dog feature, to make sure that every treatment is precise.

**One of the best  
Ultrasound technology  
for most effective  
treatment**

Frequencies : 1 Mhz, 3 Mhz & Dynamic from single applicator.

Clinical library: 25 programs.

Treatment Head : 5 cm<sup>2</sup>

Back Light Graphic LCD Display.

Duty cycles : 10%, 20%, 50% & Continuous.

Pulse Rate : 16, 48, 100 Hz

Intensity : 2.5 w/cm<sup>2</sup> Continuous and 3 w/cm<sup>2</sup> Pulse



Timer: Up to 99 Min.